



### Disposable Endoscopic Hemoclip




### Disposable Injection Needle




### Disposable Polypectomy Snare



### Disposable Biopsy Forceps



Rotatable



#### Key Features

- 360°rotatable clips design, offering a precise placement
- Atraumatic tip prevents the endoscope from damage
- Sensitive release system, easy to release clips provision
- Repeated opening and closing clip for accurate positioning
- Ergonomically shaped handle; User friendly

#### Key Features

- Strong PTFE sheath construction facilitates advancement through difficult pathways
- Needle tip angel 30°- sharp puncture
- Transparent inner tube- Can be used to observe blood return
- Shield needles designed to prevent damage to working channel
- Ergonomic handle design, easy to control the needle moving

#### Key Features

- Wire in a braided construction makes the polyps not easy to slip off
- Smooth open and close mechanism for optimum ease of use
- Made by rigid medical stainless-steel offering a precise & quick cutting properties
- Smooth sheath to prevent damage to your endoscopic channel
- Standard power connection, compatible with all main high-frequency devices on the market

#### Key Features

- Distinct catheter and position markers for visibility during insertion and withdrawal
- Coated with super-lubricious PE for better glide and protection for endoscopic channel
- Medical stainless steel, four-bar-type structure makes sampling safer and more effective
- Ergonomic handle, easy to operate
- Spike type is recommended for soft sliding tissue sampling

Model	Clip opening size(mm)	Working Length(mm)	Endoscopic Channel(mm)	Characteristics
ZRH-HCA-165-9-L	9	1650	2.8	Gastroscopy
ZRH-HCA-165-12-L	12	1650	2.8	
ZRH-HCA-165-15-L	15	1650	2.8	
ZRH-HCA-235-9-L	9	2350	2.8	
ZRH-HCA-235-12-L	12	2350	2.8	Colonoscopy
ZRH-HCA-235-15-L	15	2350	2.8	
ZRH-HCA-165-9-S	9	1650	2.8	Gastroscopy
ZRH-HCA-165-12-S	12	1650	2.8	
ZRH-HCA-165-15-S	15	1650	2.8	Colonoscopy
ZRH-HCA-235-9-S	9	2350	2.8	
ZRH-HCA-235-12-S	12	2350	2.8	
ZRH-HCA-235-15-S	15	2350	2.8	

Type	Sheath O.D. Dia. I (mm)	Working Length Lx 50 (mm)	Needle Size (diameter / length)	Endoscopic Channel (mm)
ZRH-PN-2418-214	Φ2.4	1800	21G, 4mm	2.8
ZRH-PN-2418-234	Φ2.4	1800	23G, 4mm	2.8
ZRH-PN-2418-254	Φ2.4	1800	25G, 4mm	2.8
ZRH-PN-2418-216	Φ2.4	1800	21G, 6mm	2.8
ZRH-PN-2418-236	Φ2.4	1800	23G, 6mm	2.8
ZRH-PN-2418-256	Φ2.4	1800	25G, 6mm	2.8
ZRH-PN-2423-214	Φ2.4	2300	21G, 4mm	2.8
ZRH-PN-2423-234	Φ2.4	2300	23G, 4mm	2.8
ZRH-PN-2423-254	Φ2.4	2300	25G, 4mm	2.8
ZRH-PN-2423-216	Φ2.4	2300	21G, 6mm	2.8
ZRH-PN-2423-236	Φ2.4	2300	23G, 6mm	2.8
ZRH-PN-2423-256	Φ2.4	2300	25G, 6mm	2.8

Type	Loop Width D x 20% (mm)	Working Length L x 10% (mm)	Sheath O.D. Dia. I (mm)
Oval	15/25	1200	Φ1.8
	15/25/35	1800	Φ2.4
	15/25/35	2300	Φ2.4
Hexagonal	15/25	1200	Φ1.8
	15/25/35	1800	Φ2.4
	15/25/35	2300	Φ2.4
Crescent	15/25	1200	Φ1.8
	15/25	1800	Φ2.4
	15/25	2300	Φ2.4

Model	Jaw open size (mm)	O.D. (mm)	Length (mm)	Serrated Jaw	Spike	PE Coating
ZRH-BFA-2416-PWS	6	2.4	1600	NO	NO	YES
ZRH-BFA-2423-PWS	6	2.4	2300	NO	NO	YES
ZRH-BFA-1816-PWS	5	1.8	1600	NO	NO	YES
ZRH-BFA-1812-PWS	5	1.8	1200	NO	NO	YES
ZRH-BFA-1806-PWS	5	1.8	600	NO	NO	YES
ZRH-BFA-2416-PZS	6	2.4	1600	NO	YES	YES
ZRH-BFA-2423-PZS	6	2.4	2300	NO	YES	YES
ZRH-BFA-2416-CWS	6	2.4	1600	YES	NO	YES
ZRH-BFA-2423-CWS	6	2.4	2300	YES	NO	YES
ZRH-BFA-2416-CZS	6	2.4	1600	YES	YES	YES
ZRH-BFA-2423-CZS	6	2.4	2300	YES	YES	YES

